Phenylalanine Reagent Droppers

Reagent for the phenylalanine deamination test,
in dropper bottles.

INTENDED PURPOSE

Reagent used for the differentiation of microorganisms on the basis of the capability to produce phenylpyruvic acid from phenylalanine deamination. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Phenylalanine Reagent is a solution of ferric chloride used for microbic differentiation and identification on the basis of the capability of microorganisms to produce phenylpyruvic acid from phenylalanine deamination.

KIT CONTENT

- 20 dropper bottles of Phenylalanine Reagent Droppers (2.5 ml)

METHOD PRINCIPLE

Addition of Phenylalanine Reagent Droppers to a microbial culture determines the reaction of any phenylpyruvic acid with the ferric salt to develop a green colour. The reaction is detected visually and is considered a positive test. Most members of the genera *Proteus* and *Providencia* are positive whereas most *Escherichia* and *Citrobacter* species are negative.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, test tube, culture media, physiological solution (0.85% saline), Phenylalanine Agar, quality control organisms.

REAGENTS

- Each bottles contain 10%(w/v) ferric chloride in aqueous solution.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory.

Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Inoculate the microorganism onto a slant of Phenylalanine Agar.
2. Incubate at $35 \pm 2^\circ\text{C}$ for 24 h.
3. Add 4-5 drops of Phenylalanine Reagent Droppers directly to the growth on the slant.
4. The presence of phenylpyruvic acid is indicated by the development of a green colour diffusing into the agar within 5 minutes.

INTERPRETING RESULTS

The test is positive if a visible green color develops upon addition of a solution of Phenylalanine Reagent Droppers. In case of negative reaction colour not change.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL

Appearance of Phenylalanine Reagent Droppers: Dark yellow solution.

Control strain		Incubation	Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	Up to 24 h/ 35°C ± 2 °C	Negative reaction: No color change
<i>Klebsiella pneumoniae</i>	WDCM 0097 ATCC® 13883 NCTC 9633		
<i>Staphylococcus aureus</i>	WDCM 00034 ATCC® 25923 NCTC 12981		
<i>Pseudomonas aeruginosa</i>	WDCM 00025 ATCC® 27853 NCTC 12903		
<i>Proteus mirabilis</i>	ATCC® 25933		Positive reaction: Development of a green colour
<i>Morganella morganii</i>	WDCM 25830 ATCC® 25830 NCTC 235		

PERFORMANCE CHARACTERISTICS

Performance testing of Phenylalanine Reagent Droppers was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

A heavy growth of organisms on the slant is essential. The reaction must be read promptly.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Phenylalanine Reagent Droppers	20 x 1 ml	87004

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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