## Easy Dry™ Plate Count + TTC

Medium for the enumeration of bacteria in food, water and other materials, according to APHA and ISO 4833.

### DESCRIPTION

Liofilchem Easy Dry™ are absorbent pads impregnated with a sterile, dehydrated culture medium. Each pad is preplated in a Petri dish and is immediately ready to use after pouring sterile distilled or deionized water on it. Easy Dry™ are optimal for the examination of large sample volumes by the membrane filter method.

Easy Dry™ Plate Count + TTC is a medium used for the determination of the total microbial content from food and animal feed, water and other materials.

This medium is based on the formula recommended by the American Public Health Association and ISO 4833, modified by the addition of TTC.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Casein	5.0
Yeast Extract	2.5
Glucose	1.0
Triphenyltetrazolium Chloride (TTC)	0.1
Final pH 7.0 ± 0.2 at 25°C	

### METHOD PRINCIPLE

Enzymatic digest of casein provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Glucose is the fermentable carbohydrate. 2,3,5-triphenyl tetrazolium chloride (TTC) is a redox dye which is reduced to the insoluble red formazan by actively growing cells.

### PREPARATION

1. Cut open a bag and remove the number of Easy Dry™ plates needed.
2. Moisten the pad contained in the Petri dish with 2.2 ml of sterile distilled or deionized water.
3. Wait 5 minutes before using.

### TEST PROCEDURE

Filter the sample trough a filter membrane (0.45 µm pore diameter). Transfer the membrane onto a plate containing a just rehydrated pad. Incubation conditions may vary depending on the organisms under study. For a general aerobic count, incubate aerobically at 30°C for 72 hours.

### INTERPRETING RESULTS

Examine the membrane for growth. Most colonies will be red due to reduction of TTC. Count all colonies, regardless of pigmentation, to obtain the total aerobic count. Report the count as CFU per ml of sample allowing for dilution factors.

### APPEARANCE OF THE MEDIUM

Whitish pad. Yellowish once rehydrated.

### STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

### SHELF LIFE

2 years.

## QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Incubation conditions: aerobically at  $30 \pm 1^\circ\text{C}$  for  $72 \pm 3$  hours.

### QC Table.

Microorganism		Growth	Specification
<i>Bacillus subtilis</i>	WDCM 00003	Good	Red colonies
<i>Enterococcus faecalis</i>	WDCM 00009	Good	Red colonies
<i>Escherichia coli</i>	WDCM 00012	Good	Red colonies
<i>Staphylococcus aureus</i>	WDCM 00034	Good	Red colonies
<i>Pseudomonas aeruginosa</i>	WDCM 00024	Good	Red colonies

## WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

## DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

## BIBLIOGRAPHY

- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 4833 (2003) Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of microorganisms – Colony count technique at  $30^\circ\text{C}$ .
- Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17<sup>th</sup> ed. American Public Health Association, Washington, D.C.
- Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4<sup>th</sup> ed. American Public Health Association, Washington D.C.
- Greenberg A.E, L.S. Clesceri and A.D. Eaton (1992) Standards methods for the examination of water and wastewater, 18<sup>th</sup> ed. American Public Health Association, Washington D.C.

PRESENTATION	Contents	Ref.
Easy Dry™ Plate Count + TTC    60 mm ready-to-use plates	100 plates	87501

## TABLE OF SYMBOLS

<b>LOT</b> Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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- EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 4833 (2003) Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of microorganisms – Colony count technique at  $30^\circ\text{C}$ .
- Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17th ed. American Public Health Association, Washington, D.C.
- Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington D.C.
- Greenberg A.E, L.S. Clesceri and A.D. Eaton (1992) Standards methods for the examination of water and wastewater, 18th ed. American Public Health Association, Washington D.C.

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