Nitrate Test

Rapid test for determining the ability of an organism to reduce nitrate to nitrite or free nitrogen.

INTENDED PURPOSE

Rapid test for the evaluation of reduction from nitrates to nitrites on bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Nitrate Test is used to determine the ability of an organism to reduce nitrates to nitrites. This test represents a useful component of biochemical test batteries for identification of bacteria, including differentiating members of the family Enterobacteriaceae from other gram-negative bacilli, identifying species of *Neisseria* and separating them from *Moraxella* and *Kingella* species, as well as facilitating species identification of *Corynebacterium* and other non-sporing gram-positive rods.

KIT CONTENTS

- 30 Tubes with a dehydrated culture medium containing Potassium Nitrate
- 1 Bottle of N,N-Dimethyl-alpha-naphthylamine reagent (3 ml)
- 1 Bottle of Sulphanilic acid reagent (3 ml)

METHOD PRINCIPLE

Reduction of nitrate (NO_3) to nitrite (NO_2), which is mediated by the bacterial enzyme nitrate reductase, indicates that the organism can use nitrate as an electron acceptor during anaerobic respiration. Nitrites are colourless and react with the reagents sulphanilic acid and alpha-naphthylamine to form a red compound.

Positive reaction: the microorganisms reduce nitrates to nitrites resulting in a red-orange color development after the addition of reagents to the culture medium.

Negative reaction: microorganisms that are unable to reduce nitrates to nitrites do not produce a colour development after the addition of the two reagents to the culture medium.

However, the formation of other end products of nitrate reduction, including nitric oxide (NO), nitrous oxide (N_2O) and ammonia (NH_3), that do not cause a color change of the medium, is also possible; Zinc dust is therefore used to detect unreduced nitrate. No color change after the addition of zinc dust to the medium indicates that nitrate has been reduced to nitrogen gas (N_2) and the organism is nitrate positive. The appearance of a red color after the addition of zinc dust to a colorless reaction tube indicates a negative reaction, i.e., the organism has failed to reduce NO_3 .

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: zinc dust, incubator, inoculating loop, pipettes, culture media, physiological solution (0.85% saline), test tubes, 0.5 McFarland turbidity standard, quality control organisms.

REAGENTS

- Dehydrated basal medium with Potassium Nitrate
- 0.6% N,N-Dimethyl-alpha-naphthylamine in 5N Acetic acid.
- 0.8% Sulphanilic acid in 5N Acetic acid.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Allow the kit components to reach room temperature before use.
2. Add 0.2 mL of physiological solution to a test tube.
3. Inoculate heavily with the test organism(*) from a fresh overnight pure culture.
4. Add 0.2 mL of sterile physiological solution to another tube. Do not inoculate this tube and use it as negative control;

5. Incubate both inoculated and uninoculated tubes at $35 \pm 2^\circ\text{C}$ for a minimum of 2 hours and up to 24 h(**).
6. Add one drop of Alpha-naphthylamine reagent and one drop of Sulphanilic acid reagent to each tube. Shake gently to mix reagents.
7. Observe for a red-orange colour development within 2 minutes.
8. In case of a negative test result (no color development) add a small quantity of zinc dust to the medium.
9. Observe for a red-orange colour development within 10 minutes.

* Microbial suspension equivalent to 0.5 McFarland turbidity standard (which corresponds to $1-2 \times 10^8$ CFU/ml for *E. coli*).

** Prolonged incubation may be required for some organisms to achieve a positive reaction.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

A positive nitrate reaction is indicated by either a colour change to orange-red after within 2 minutes addition of N,N-Dimethyl-alpha-naphthylamine and Sulphanilic acid reagents OR no orange-red colour development within 10 minutes after addition of zinc dust.

If the organism is nitrate negative there is no a color development after addition of N,N-Dimethyl-alpha-naphthylamine and Sulphanilic acid reagents AND an orange-red color appears after addition of zinc dust.

Further tests are required for confirmation.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year

QUALITY CONTROL

Expected Response:

Control strain	Characteristic reactions	
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	Positive reaction: Color change to orange-red
<i>Listeria monocytogenes</i>	WDCM 00109 ATCC® 35152 NCTC 7973	Negative reaction: No color change

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Nitrate Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

The test of nitrate reduction to nitrite is helpful for the bacterial identification. It must be followed yet by other tests as the Gram colouring, biochemical and serological tests, etc.

Interpretation of colour reactions should be made soon after addition of reagents.

Zinc dust must be added to any negative reaction to eliminate the possibility of a false negative reaction due to the organism's ability to reduce the nitrate to nitrogen gas.

The addition of an excessive quantity of zinc dust could result into false negative results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**

- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Nitrate Test	30 tests	88009

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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
References / Riferimenti

1. Ewing 1986. Edwards & Ewing's. Identification of Enterobacteriaceae, 4th ed. Elsevier Science Publishing Co., Inc., New York.
2. MacFaddin 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore
3. Finegold and Baron. 1986 Bailey and Scott's. Diagnostic microbiology, 7th ed. The C.V. Mosby Company, St. Louis.
4. Kelly, Brenner and Farmer. 1985. In Lennette, Balows, Hausler and Shadomy (ed.), Manual of clinical microbiology, 4th ed. ASM, Washington, D.C.

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