Urea/Indole Test

Rapid test for detection of urease and indole production.

INTENDED PURPOSE

Rapid test for the detection of urease and the evaluation of indole production used for the differentiation and identification of bacterial strains. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Urea/Indole Test is used for the rapid qualitative determination of urease activity and indole production as a first step in screening non-lactose fermenting colonies from agar plates.

KIT CONTENT

- 30 Tubes of Urea/Indole Test
- 1 Tube of Kovac's Reagent (3.0 ml)

METHOD PRINCIPLE

The enzyme urease catalyzes the hydrolysis of urea to carbon dioxide, water and ammonia with increase in pH and a color change of the phenol red indicator to red-fuchsia. Bacteria that produce the enzyme tryptophanase are able to degrade the amino acid tryptophan into pyruvic acid, ammonium and indole. The indole produced reacts with p-dimethylaminobenzaldehyde of the Kovac's (indole) reagent to form a red-violet compound.

Enteric pathogens such as *Salmonellae* and *Shigellae* are urease negative whilst organisms such as *Proteus*, *Morganella*, *Klebsiella* spp and some species of *Citrobacter* are strongly urease positive. The production of indole from tryptophan is a characteristic absent in *Salmonella* but present in *E. coli*, *Morganella* and some species of *Klebsiella*, *Aerobacter* and *Citrobacter*.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, culture media, physiological solution (0.85% saline), test tubes, 0.5 McFarland turbidity standard, quality control organisms.

REAGENTS

- A tube of Urea/Indole Test contains desiccated biochemical substrates and nutritive components such as urea, phenol red and tryptone, which includes tryptophane.
- Kovac's Reagent consists of 5% p-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloric acid and 75% isobutyl alcohol.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Urea/Indole Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.3 mL of physiological solution to each tube.
3. Inoculate heavily with the test organism(*) from a fresh overnight pure culture.
4. Incubate at $35 \pm 2^\circ\text{C}$ for up to 24 hours(**).
5. Examine tubes at intervals of 1, 4 and 24 hours for broth colour change.
6. Add 2-3 drops of Kovac's Reagent to perform the Indole Test.

* Microbial suspension equivalent to 0.5 McFarland turbidity standard (which corresponds to $1-2 \times 10^8$ CFU/ml for *E. coli*).

** Prolonged incubation may result in a false-positive test due to hydrolysis of proteins in the medium.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

A positive urease reaction is indicated by a colour change to red-fuchsia.

The appearance of a red ring at the surface of the broth within 30 seconds (after addition of Kovac's Reagent) indicates a positive indole reaction.

If the organism is urease negative and indole negative it is possibly *Salmonella*.

Further tests are required for confirmation.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of Urea/Indole Test: Yellowish powder, yellow solution once reconstituted.

Appearance of Kovac's Reagent: Yellow-green liquid, clear.

Control strain		Inoculum	Incubation	Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241			Negative urease reaction: No color change Positive indole reaction: Formation of a red ring
<i>Proteus mirabilis</i>	WDCM 00023 ATCC® 29906 NCTC 11938	1-2 x 10 ⁸ CFU/ml	Up to 24 h/ 35 ± 2°C	Positive urease reaction: Colour change to red fuchsia Negative indole reaction: Yellow/brown ring
<i>Salmonella</i> Typhimurium	WDCM 00031 ATCC® 14028 NCTC 12023			Negative urease and indole reactions

PERFORMANCE CHARACTERISTICS

Performance testing of Urea/Indole Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium providing false negative results.

As the urea test is based on demonstration of alkalinity, it is not specific for urease. The hydrolysis of proteins by contaminating organisms in the medium may raise the pH to alkalinity resulting in false positive reactions.

Allow at least 30 seconds for the color to develop before considering the indole test negative. Pigmented organisms may yield ambiguous results.

Indole produced by certain anaerobes, such as *Clostridium* species, can be rapidly degraded resulting in false negative reactions.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.

- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Urea/Indole Test	30 tests	88024

Revision History

Revision	Release Date	Change Summary
3	2024-09-16	Added notice to report any malfunction, defect or incident.
2	2023-04-24	Updated layout and content in compliance with IVDR 2017/746.

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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
lio.filchem.com/ifu-sds

References / Bibliografia / Bibliografia

- EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- Public Health England. (2019). Urease test. UK Standards for Microbiology Investigations. TP 36 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
- Public Health England. (2018). Indole test. UK Standards for Microbiology Investigations. TP 19 Issue 4
- Isemberg H.D. and Sundheim L.H.(1980) Indole reactions in bacteria. J. Bacteriol. 75:682-690. Baltimore, Williams & Wilkins.
- MacFaddin J.F. (1980) Biochemical tests for identification for medicalbacteria. 2nd ed. 173-183. Baltimore, Williams & Wilkins.
- Blazevic D.J. and Ederer, G.M.(1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

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