O129 Disc

Diagnostic discs for differentiation of *Vibrio* species.

INTENDED PURPOSE

Test for differentiation of *Vibrio* spp based on sensitivity to vibrio-static agent O129 (2,4-diamino-6,7-di-isopropyl-pteridine phosphate). This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

O129 Disc is a test used for differentiation of *Vibrio* spp based on sensitivity to vibrio-static agent O129 (2,4-diamino-6,7-di-isopropyl-pteridine phosphate).

KIT CONTENT

- 1 cartridge of 50 discs.

Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

METHOD PRINCIPLE

This test is useful in the differentiation of *Vibrio* from other gram-negative bacteria especially *Aeromonas*, which are characteristically resistant to O129. Even among the genus *Vibrio*, different species show different sensitivity to O129. Hence two different concentration discs can be simultaneously tested to determine the degree of sensitivity of the species.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, culture media, physiological solution (0.85% saline), 0.5 McFarland turbidity, forceps, quality control organisms.

REAGENTS

- O129 Disc is available at two concentrations: 10 µg and 150 µg

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Inoculate the plate* by streaking the sample on the agar surface to produce even growth.
2. Using sterile forceps, place one 150 µg disc only onto the agar surface to differentiate between *Vibrio* spp and other gram-negative bacteria. Otherwise use both the 10 µg disc and the 150 µg disc to also differentiate within *Vibrio* species. Ensure suitable distance between the discs to avoid zones of inhibition overlapping.
3. Incubate aerobically at 35°C for 24 hours.
4. Observe for zones of inhibition around the disc.

***Note:** It is recommended to use non-selective salt containing media (0.5%) such as Nutrient Agar or Columbia Agar (Horse Blood 5%).

INTERPRETING RESULTS

Most *Vibrio* species are sensitive with 150 µg but species differ with 10 µg discs (some strains of *V. cholerae* O1 and O139 may be resistant to both disc contents). Test organism could result sensitive (S), partially sensitive (PS) or resistant (R), to O129 as indicated in table 1.

Table 1.

Organism	Zone of inhibition around the disc		Result
	O129 10 µg	O129 150 µg	
<i>Vibrio cholerae</i>	Yes	Yes	S
<i>Vibrio vulnificus</i>	Yes	Yes	S
<i>Vibrio</i> spp	Yes	Yes	S
<i>Vibrio parahaemolyticus</i>	No	Yes	PS
<i>Vibrio alginolityticus</i>	No	Yes	PS
<i>Vibrio fluvalis</i>	No	Yes	PS
<i>Vibrio furnissii</i>	No	Yes	PS
<i>Vibrio cholerae</i> O1	No	No	R
<i>Vibrio cholerae</i> O139	No	No	R
<i>Aeromonas</i> spp	No	No	R

Sensitive (S): zone of inhibition around both O129 discs.

Partially Sensitive (PS): zone of inhibition around the 150 µg disc and no zone around 10 µg disc.

Resistant (R): no zone of inhibition around both O129 discs

STORAGE

2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Control strain		Inoculum	Incubation	Characteristic reactions
<i>Vibrio parahaemolyticus</i>	ATCC® 17802	McFarland 0.5	24 h, 35 ± 2°C	Inhibition zone
<i>Aeromonas hydrophila</i>	ATCC® 7966	McFarland 0.5	24 h, 35 ± 2°C	No inhibition zone

PERFORMANCE CHARACTERISTICS

Performance testing of O129 Disc was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Negative O129 tests must be interpreted with caution. Strains of *Vibrio cholerae* resistant to O129 and trimethoprim have been reported. Both plasmids and transposons appear to be responsible. Specialised sensitivity testing media which contain low levels of sodium chloride should not be used as this may give poor growth of halophilic vibrios and may also give false positive results with Enterobacteriaceae. Nutrient Agar can be recommended for such testing. It is recommended that the tests are performed on colonies from pure culture.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.

- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Concentration	Packaging	Ref.
O129 Disc	150	50 discs	88043
O129 Disc	10	50 discs	88044

Revision History

Revision	Release Date	Change Summary
2	2024-09-16	Added notice to report any malfunction, defect or incident.
1	2023-05-15	Updated layout and content in compliance with IVDR 2017/746.

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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