DESCRIZIONE

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F01013
Rev.4 / 04.11.2016

DESCRIPTION

Dermatest is a ready-for-use system for dermatophyte fungus detection, with culture media solidified onto a plastic support. Such systems are universally known as dip-slide. The results can be used as an aid in diagnosis of skin infections and must be used in conjunction with data on the patient's clinical history.

CONTENTS OF THE PACKAGES

Dermatest is available in the packages of 10 and 120 tests.

CONFIGURATION

Dermatest

- | | | |
|-----------------------------|----------------------------------|--|
| 1. Malt Agar (amber) | 2. Corn Meal Agar (white) | 3. Dermatophyte Test Medium Agar (orange) |
|-----------------------------|----------------------------------|--|

PRINCIPLE OF THE METHOD

Dermatest inoculated with a collection swab containing clinical material permits a search for dermatophyte fungi. Dermatophyte Test Medium Agar is prepared according to formula of Taplin, Zaias and Rebell and permits rapid selective isolation of dermatophyte fungi responsible for lesions of the skin, nails and hair such as *Microsporum* spp., *Trichophyton* spp. and *Epidermophyton* spp. The selection agents in the medium, cycloheximide, gentamycin and chloramphenicol, totally eliminate bacterial contamination and markedly reduce contamination by saprophytic fungi, at the same time permitting a high rate of isolation of dermatophytes. Identification of dermatophytes on these media is based on the use of phenol red as pH indicator. The dermatophytes produce alkaline metabolites which turn the indicator from yellow to red. The reliability of this characteristic as a differentiating mechanism has been confirmed in various studies by Allen, Carrol, Taplin *et al.*, who report an accuracy of 97-100% in identification of dermatophytes with Dermatophyte Test Medium Agar. The medium makes it possible to perform a dermatophyte diagnosis following at least 48 hours of incubation. Malt Agar is used to isolate opportunistic pathogenic fungi such as *Aspergillus* spp., *Fusarium* spp., *Mucor* spp. and *Rhizopus* spp., from samples of clinical origin for the isolation of pathogens sensitive to cycloheximide (*Allescheria boydii* and *Cryptococcus neoformans*). Corn Meal Agar is a suitable medium for the isolation of *Candida albicans*. When colonies of *Candida albicans* grown on this medium are observed under the microscope, one notes the presence of characteristic chlamydospores, the accepted criterion for morphological identification of this species (Prospero and Reyes).

TYPICAL FORMULATIONS FOR THE MEDIA USED ON DERMATEST (grams/litre)

MALT AGAR		CORN MEAL AGAR		DTM AGAR	
Malt extract	30.0	Corn meal extract	2.0	Soya peptone	10.0
Agar	17.0	Agar	15.0	Glucose	10.0
pH 5.5 ± 0.2		pH 6.0 ± 0.2		Phenol red	0.2
				Cycloheximide	0.5
				Gentamicin	0.1
				Chloramphenicol	0.1
				Agar	20.0
				pH 5.5 ± 0.2	

COLLECTION AND CONSERVATION OF SAMPLES

The material taken from lesions of the skin, hair and nails must be seeded onto the surface of the two media, ensuring good adherence to the media.

TEST PROCEDURE

1. Take a dip-slide from the container
2. Identify the sample with an adhesive label or by writing the data directly on the dip-slide container.
3. Unscrew and extract the slide from the cylindrical container without touching the surface of the culture media.
4. Holding the slide by the stopper, seed the surfaces of the culture media directly with a collection swab containing clinical material taken from the cutaneous lesions etc.
5. Put the slide back into the container and screw it shut.
6. Incubate aerobically at 25-30°C for 48 hours for up to a maximum of 7 days.

INTERPRETATION OF THE RESULTS

Following incubation, the presence of dermatophytes is revealed by the appearance of colonies on the surface of the culture media. On Dermatophyte Test Medium Agar the dermatophyte fungi develop with red colouring around the colonies. In the case of fast growing dermatophytes, the red colour appears after 48 hours of incubation. In the case of slow growing dermatophytes, from 3 to 7 days of incubation are necessary. If the colonies are small, the red colour is restricted to the immediate vicinity of the colonies; when they merge and are conspicuous, the indicator changes colour

over the whole medium. Saprophytic fungi do not normally grow on Dermatophyte Test Medium Agar unless the sample is very heavily contaminated; in such a case, growth is later, with characteristic colours for the colonies (black in the case of *Aspergillus niger* and *Cladosporium* spp., green in the case of *Penicillium* spp.) at times with the medium turning red: *Candida albicans* develops without causing any change in the colour of the medium. On Malt Agar and on Corn Meal Agar, the same micro-organisms grow with the colours characteristic of the filamentous fungus isolated.

QUALITY CONTROL

For quality control of the Dermatest, the dip-slides are inoculated with a broth culture of microbes listed in the table. Conditions for incubation: 25-30°C for 48-96 h under aerobic atmosphere.

Strains used and expected results in quality control:

Microorganism	MALT AGAR	CORN MEAL AGAR	DTM AGAR
<i>Escherichia coli</i> ATCC® 25922	colourless colonies	colourless colonies	growth inhibited
<i>Pseudomonas aeruginosa</i> ATCC® 27853	colourless colonies	colourless colonies	growth inhibited
<i>Candida albicans</i> ATCC® 10231	white colonies	white colonies	white colonies (no color change in the medium)
<i>Trichophyton mentagrophytes</i> ATCC® 9533	white colonies	white colonies	white colonies (red medium)
<i>Aspergillus niger</i> ATCC® 16404	black colonies	black colonies	growth inhibited

PRECAUTIONS

The media described here are not classified as hazardous by current legislation, nor do they contain harmful substances in concentrations $\geq 1\%$. No safety data sheets are therefore required for them. Dermatest are disposable devices.

Dermatest is only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in a laboratory by adequately trained personnel using approved asepsis and safety methods for dealing with pathogenic agents.

DIP-SLIDE CONSERVATION

Store at 10-25°C in the original packaging, which provides protection against light. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions, the products described here will be valid until the expiry date shown on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration (conspicuous changes in the colour of the media, drying or detachment of the medium from its support, presence of fungal or bacterial colonies). DO NOT FREEZE.

DISPOSING OF USED MATERIAL

After use, Dermatest and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

1. Prospero Magdalene T. and Reyes A.C. (1955) Acta Mel. Phillipina 12 (2). 69-74.
2. Allen A.M., Drewry R.A., Weaver R.E. (1970) Arch. Derm. 102, 68-70.
3. Carrol H.F. (1974) J.A.V.M.A. 155, 192-195.
4. Taplin D., Zaias N., Rebell G., Blank H. (1969) Arch. Derm., 99, 1969.

PRESENTATION

Dermatest	REF	50021	51021
		10 Dip Slide	120 Dip Slide

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM® S.r.l.