Tryptose Phosphate Broth

Liquid medium for cultivating fastidious organisms from clinical and nonclinical specimens.

TYPICAL FORMULA	(g/l)
Enzymatic Digest of Casein	20.0
Glucose	2.0
Sodium Chloride	5.0
Disodium Phosphate	2.5
Final pH 7.3 ± 0.2 at 25°C	

DESCRIPTION

Tryptose Phosphate Broth is a liquid medium used for the cultivation of a wide variety of fastidious organisms from clinical specimens and other materials.

PRINCIPLE

Enzymatic digest of casein provides amino acids, nitrogen, carbon, minerals, vitamins and minerals for the growth of microorganisms. Glucose is a source of energy. Sodium chloride maintains the osmotic balance of the medium. Disodium phosphate is the buffering agent.

PREPARATION

Suspend 29.5 g of powder in 1 liter of deionized or distilled water (*). Bring to boil and shake until completely dissolved. Dispense in final containers. Sterilize at 121°C for 15 minutes.

(*) Add 1-2 g of agar to facilitate anaerobic growth, especially in the lower zone of the medium.

TECHNIQUE

Inoculate the tube with the specimen. Incubate at 35 ± 2°C for 18-48 hours.

INTERPRETATION OF RESULTS

Turbidity indicates microbial growth.

STORAGE AND TRANSPORT CONDITIONS

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store the prepared medium at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- Harmon, S. M., D. A. Kautter, D. A. Golden and E. J. Rhodehamel (1995) FDA Bacteriological analytical manual, 8th ed. AOAC International, Arlington, VA.
- MacFaddin J. D. (1985) Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1, p. 802-804. Williams & Wilkins, Baltimore, MD.
- Ginsberg, Gold and Jordan (1955) Proc Soc Exp Biol Med. 89:66.

PRODUCT SPECIFICATIONS

NAME

Tryptose Phosphate Broth

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
610364	500 g	500 g of powder in plastic bottle
620364	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.3 ± 0.2

USE

Tryptose Phosphate Broth is a liquid medium used for the cultivation of a wide variety of fastidious organisms from clinical specimens and other materials of sanitary importance

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

Ready-to-use medium

Appearance: clear to very slightly opalescent

Colour: Light amber

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: ≤100 CFU
Incubation Conditions: 18-48 h at 35 ± 2°C

Microorganism

Microorganism	ATCC®	Growth
<i>Escherichia coli</i>	ATCC® 25922	Good
<i>Neisseria meningitidis</i>	ATCC® 13090	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Good
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Good
<i>Streptococcus pyogenes</i>	ATCC® 19615	Good

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.