## COAGULASE TEST

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## COAGULASE TEST

Freeze-dried rabbit plasma for coagulase test.

### TYPICAL FORMULA

Freeze-dried rabbit plasma with EDTA..... 4.0 ml

### DESCRIPTION

COAGULASE TEST is constituted by rabbit freeze-dried plasma with EDTA, for the determination of enzyme coagulase produced by *Staphylococcus aureus*.

### PRINCIPLE

The test shows the coagulant activity of potentially pathogenic staphylococci on plasma. This activity is due to two factors at least: free coagulase ( extracellular enzyme) and bound coagulase or clumping factor (antigen of cellular wall).

The enzyme is present in the most of biotypes belonging to the species of *Staphylococcus aureus* and in biotypes belonging to the species *S. intermedius* and *S. hyicus*, pathogenic opportunist for animals, and is always absent in saprophyte and commensal species.

### TECHNIQUE

1. Take a bottle of freeze-dried plasma from the kit and aseptically reconstitute with 4 ml of sterile physiological solution. Mix until completely dissolved, avoiding foam formation.
2. Prepare a culture broth in Brain Heart Infusion, taking one or more colonies from selective culture media for *Staphylococcus aureus* isolation, and incubate at 36+/-1 °C for 4-6 hours.
3. Mix 0.5 ml of Coagulase Test with 0.5 ml of culture broth and incubate at 36+/-1 °C for 1-2-4-8-24 hours.
4. Examine for the formation of coagulum, carefully inclining the tube to one side, without shaking. Do not incubate after 24 hours because cases of fibrinolysis may occur.

### INTERPRETATION OF RESULTS

In the most or in the whole culture medium the presence of enzyme coagulase is shown by the formation of a well defined coagulum. The test determines both free and bound coagulase.

### STORAGE

2-8 °C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING and PRECAUTIONS

The product contains dangerous substances according to directives 1999/45/CE and 2001/60/CE or for which exist recognized exposure limits.

### DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.

### REFERENCES

1. W.E. Kloos and J.H. Jorgensen "Staphylococci" p. 143-153. In E.H. Lennette, A. Balows, W.J. Hausler, H.J. Shadomy (eds) *Manual of Clinical Microbiology*, 4<sup>th</sup> Edition, American Society for Microbiology, Washington, D.C. 1985.

## PRODUCT SPECIFICATIONS

### NAME

COAGULASE TEST

### PRESENTATION

Glass bottles containing 4 ml rabbit freeze-dried plasma with EDTA.

### STORAGE

2-8 °C

### PACKAGING

| Code  | Content             | Packaging          |
|-------|---------------------|--------------------|
| 88030 | 5 bottles x 40 test | • 5 bottles in box |

### USE

COAGULASE TEST is constituted by rabbit freeze-dried plasma with EDTA, for the determination of enzyme coagulase produced by *Staphylococcus aureus*.

### TECHNIQUE

Refer to technical sheet of the product.

### APPEARANCE of the MEDIUM

Compact freeze-dried, beige to light pink.

### SHELF LIFE

2 years

### QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control  
Incubation conditions: 12-24 h at  $36 \pm 1$  °C, in aerobiosis

| Microorganisms               |            | Coagulation |
|------------------------------|------------|-------------|
| <i>Escherichia coli</i>      | ATCC 25922 | -           |
| <i>Staphylococcus aureus</i> | ATCC 25923 | +           |

### TABLE OF SYMBOLS

| Symbol | Meanings                           |
|--------|------------------------------------|
|        | Catalogue number                   |
|        | In vitro Diagnostic Medical Device |
|        | Manufacturer                       |
|        | Temperature limitation             |
|        | Use by                             |
|        | Batch code                         |
|        | Consult accompanying documents     |