



TEST RESULTS FORM

F10123
Rev 0 / 01.12.05

HACCP SYSTEM Plus

System for determination of total microbial count and search for and presumptive identification of pathogenic microorganisms that can be isolated from work surfaces.

Surface data
Date of sample collection
Notes

Mark the obtained result in the relevant square

well	TOTAL MICROBIAL COUNT	Well colour			
		Positive reaction		Negative reaction	
1-GR+	Low microbial contamination (4-50 CFU/cm ²)	yellow	<input type="checkbox"/>	blue	<input type="checkbox"/>
2-GR++	Moderate microbial contamination (50-500 CFU/cm ²)	yellow	<input type="checkbox"/>	blue	<input type="checkbox"/>
3-GR+++	High microbial contamination (>500 CFU/cm ²)	yellow	<input type="checkbox"/>	blue	<input type="checkbox"/>

well	BIOCHEMICAL REACTIONS FOR IDENTIFICATION	Well colour			
		Positive reaction		Negative reaction	
4-LDC	Decarboxylation of lysine	red	<input type="checkbox"/>	yellow	<input type="checkbox"/>
5-H ₂ S	Production of hydrogen sulphide	black	<input type="checkbox"/>	yellow	<input type="checkbox"/>
Presence of Salmonella : 1-LDC red 2-H ₂ S black		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Presence of Citrobacter : 1-LDC yellow 2-H ₂ S black		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Presence of Proteus : 1-LDC yellow 2-H ₂ S black		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

well	PRESUMPTIVE IDENTIFICATION	Well colour			
		Positive reaction		Negative reaction	
6-PRO	<i>Proteus spp.</i> / <i>Providencia spp.</i>	brown-black	<input type="checkbox"/>	yellow	<input type="checkbox"/>
7-ESC	<i>Escherichia coli</i>	blue	<input type="checkbox"/>	grey-red	<input type="checkbox"/>
8-IND	Indole test	pink-red ring	<input type="checkbox"/>	yellow	<input type="checkbox"/>
9-PSE	<i>Pseudomonas spp.</i>	turbid green	<input type="checkbox"/>	colourless	<input type="checkbox"/>
10-LIS	<i>Listeria spp.</i>	black	<input type="checkbox"/>	yellow	<input type="checkbox"/>
11-STA	<i>Staphylococcus aureus</i>	black	<input type="checkbox"/>	yellow	<input type="checkbox"/>
12-Y/M	Yeasts and moulds	yellow	<input type="checkbox"/>	green	<input type="checkbox"/>

Date of Test	Operator
--------------	----------