


Disc kits for Antimicrobial Resistance Testing

KPC&MBL&OXA-48 discs kit


Description		µg	Packaging	Ref.
Meropenem	MRP	10	50 Test	99007
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			
Temocillin	TMO	30		



Klebsiella pneumoniae
clinical isolate
OXA-48 positive

KPC&MBL discs kit (according to EUCAST)

Description		µg	Packaging	Ref.
Meropenem	MRP	10	50 Test	99003
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			



Klebsiella pneumoniae
ATCC® BAA-1705™
KPC positive

ESBL discs kit (according to EUCAST)


Description		µg	Packaging	Ref.
Cefotaxime	CTX	30	50 Test	99004
Cefotaxime + Clavulanic acid	CTL	40 (30+10)		
Ceftazidime	CAZ	30		
Ceftazidime + Clavulanic acid	CAL	40 (30+10)		
Cefepime	FEP	30		
Cefepime + Clavulanic acid	FEL	40 (30+10)		



Disc kits for Antimicrobial Resistance Testing

ESBL discs kit (according to CLSI)


Description		µg	Packaging	Ref.
Cefotaxime	CTX	30	50 Test	99005
Cefotaxime + Clavulanic acid	CTL	40 (30+10)		
Ceftazidime	CAZ	30		
Ceftazidime + Clavulanic acid	CAL	40 (30+10)		



Klebsiella pneumoniae subsp. *pneumoniae*
ATCC® 700603™
ESBL positive

ESBL (chromosomal ind. AmpC) discs kit (according to EUCAST)

Description		µg	Packaging	Ref.
Cefepime	FEP	30	100 Test	99006
Cefepime + Clavulanic acid	FEL	40 (30+10)		



Enterobacter cloacae
clinical isolate
ESBL positive

AmpC discs kit

Description		µg	Packaging	Ref.
Cefotaxime		30	50 Test	99009
Cefotaxime + Cloxacillin				
Ceftazidime		30		
Ceftazidime + Cloxacillin				

ESBL+AmpC screen discs kit

Description		µg	Packaging	Ref.
Cefotaxime		30	50 Test	99008
Cefotaxime + Clavulanic acid		40 (30+10)		
Cefotaxime + Cloxacillin				
Cefotaxime + Clavulanic acid + Cloxacillin				

