

# Ready BAGS

Ready to use liquid media in bags of 3 or 5 liter  
for pre-enrichment of *Salmonella spp.* and enrichment of  
*Listeria monocytogenes* in foods



<p><b>LISTERIA DEMI FRASER BROTH - BAG</b> enrichment of <i>Listeria monocytogenes</i> in foods (ISO 11290)</p>	<p>499020 499025</p>	<p>3X3 L 3X5 L</p>
<p><b>LISTERIA UVM MOD BROTH - BAG</b> selective enrichment of <i>Listeria monocytogenes</i> in foods (IDF formulation)</p>	<p>499010 499015</p>	<p>3X3 L 3X5 L</p>
<p><b>BUFFERED PEPTONE WATER - BAG</b> pre-enrichment of <i>Salmonella spp.</i> in foods (ISO 6579)</p>	<p>499030 499035</p>	<p>3X3 L 3X5 L</p>
<p><b>MAXIMUN RECOVERY BROTH - BAG</b> isotonic diluent for optimal recovery of microorganisms (ISO 6887)</p>	<p>499040 499045</p>	<p>3X3 L 3X5 L</p>

 **Liofilchem® s.r.l.**

Via Scozia, Zona Industriale  
64026 Roseto degli Abruzzi, Italy  
Ph: +39 085 8930745 Fx: +39 085 8930330

[www.liofilchem.net](http://www.liofilchem.net)

F51968 Rev. 0/08.03.2012