

Mini Incubator

Description	Packaging	Ref.
MINI INCUBATOR EM4 ATL	1 unit	96750

Incubators

Description	Packaging	Ref.
M91 MICRO INCUBATOR 6 L Range of temperature: 5 to 60 °C. Capacity: 6 liters.	1 unit	96751



DIGITAL THERMOBLOCK Incubator for ampoules and tubes. Range of temperature: 5 to 150 °C.	1 unit	96752
--	--------	-------



ALUMINUM BLOCK FOR STERILTEST Aluminum block with 20 holes.	1 unit	96753
ALUMINUM BLOCK FOR STERIL CONTROL Aluminum block with 12 holes.	1 unit	96754
ALUMINUM BLOCK FOR STRIP CONTROL MiRA Test ,MeRA Test Aluminum block with 12 holes.	1 unit	96755

Air sampler



Description	Packaging	Ref.
PORTABLE AIR SAMPLER "55" Instrument for microbiological monitoring of environment with 55 mm Contact Plates, complete of: carrying case, air sampler, spare aluminum sampling head, remote control, battery charger, pencil torch for head sterilization, calibration certificate.	1 unit	96760
PORTABLE AIR SAMPLER "90" Instrument for microbiological monitoring of environment with 90 mm Petri dishes, complete of: carrying case, air sampler, spare aluminum sampling head, remote control, battery charger, pencil torch for head sterilization, calibration certificate.	1 unit	96770