



**Multidisc with eight antimicrobial agents**



Description		µg	CLSI	EUCAST	Packaging	Ref.
<b>ANAEROBES</b>	Carbenicillin	CAR	100		100 Multodiscs	95200
	Cefoxitin	FOX	30			
	Clindamycin	CD	2			
	Chloramphenicol	C	30			
	Erythromycin	E	15			
	Metronidazole	MTZ	5			
	Ampicillin	AMP	10			
	Tetracycline	TE	30			
<b>ENTEROCOCCI</b>	Tigecycline	TGC	15	✓	100 Multodiscs	95210
	Ampicillin	AMP	2	✓		
	Linezolid	LNZ	10	✓		
	Nitrofurantoin	F	100	✓		
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Vancomycin	VA	5	✓		
	Gentamicin	CN	30	✓		
	Teicoplanin	TEC	30	✓		
<b>ENTEROBACTERIA URINE</b>	Levofloxacin	LEV	5	✓	100 Multodiscs	95230
	Ampicillin	AMP	10			
	Cefuroxime	CXM	30	✓		
	Ceftazidime	CAZ	10			
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Gentamicin	CN	10	✓		
	Nitrofurantoin	F	100	✓		
	Norfloxacine	NOR	10	✓		
<b>ENTEROBACTERIA 1</b>	Amikacin	AK	30	✓	100 Multodiscs	95220
	Aztreonam	ATM	30	✓		
	Piperacillin + Tazobactam	TZP	36	✓		
	Chloramphenicol	C	30	✓		
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Netilmicin	NET	10	✓		
	Cefotaxime	CTX	5	✓		
	Nitrofurantoin	F	100	✓		
<b>ENTEROBACTERIA 2</b>	Levofloxacin	LEV	5	✓	100 Multodiscs	95240
	Cefoxitin	FOX	30	✓		
	Piperacillin + Tazobactam	TZP	36	✓		
	Cefuroxime	CXM	30	✓		
	Streptomycin	S	300			
	Netilmicin	NET	10	✓		
	Ampicillin + Sulbactam	AMS	20	✓		
	Tobramycin	TOB	10	✓		
<b>PSEUDOMONAS</b>	Amikacin	AK	30	✓	100 Multodiscs	95250
	Ticarcillin	TC	75	✓		
	Piperacillin + Tazobactam	TZP	36	✓		
	Ceftazidime	CAZ	10	✓		
	Gentamicin	CN	10	✓		
	Netilmicin	NET	10	✓		
	Levofloxacin	LEV	5	✓		
	Tobramycin	TOB	10	✓		
<b>STAPH</b>	Cefoxitin	FOX	30	✓	100 Multodiscs	95260
	Levofloxacin	LEV	5	✓		
	Netilmicin	NET	10	✓		
	Tetracycline	TE	30	✓		
	Erythromycin	ERY	15	✓		
	Gentamicin	CN	10	✓		
	Rifampicin	RD	5	✓		
	Linezolid	LNZ	10	✓		

**Multodisc with eight antimicrobial agents**



Description		µg	CLSI	EUCAST	Packaging	Ref.
<b>ACINETOBACTER</b>	Doripenem	DOR	10	✓	100 Multodiscs	95270
	Meropenem	MRP	10	✓		
	Levofloxacin	LEV	5	✓		
	Gentamicin	CN	10	✓		
	Trimethoprim-Sulfamethoxazole	SXT	25	✓		
	Tobramycin	TOB	10	✓		
	Ciprofloxacin	CIP	5	✓		
	Netilmicin	NET	10	✓		
<b>YEASTS</b>	Caspofungin	CAS	5	✓	100 Multodiscs	95280
	Fluconazole	FLU	25	✓		
	Posaconazole	POS	5	✓		
	Voriconazole	VO	1	✓		
	Amphotericin B	AMB	20			
	Ketoconazole	KCA	10			
	Flucytosine	AFY	1			
	Nystatin	NY	100 IU			