

Spore suspensions for hydrogen peroxide sterilization		
Description	Packaging	Ref.
SUSPENSION SPORE GST VHP E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82030
SUSPENSION SPORE GST VHP E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82031
SUSPENSION SPORE GST VHP E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82032
SUSPENSION SPORE GST VHP E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82033