



STERILE SWABS WITHOUT ENVELOPE

Swabs for collection and transport of samples.

DESCRIPTION

STERILE SWABS WITHOUT ENVELOPE are constituted by sterile dry swabs inserted in a tube and endowed with a cap which hermetically closes the neck of the tube, so avoiding contamination of the tube and protecting user's hands. Sterilization is assured by irradiation at 2.5 Mrad rays.

TECHNIQUE

Take a swab from the tube and use it to collect the sample. Then insert again the swab into the tube and screw the cap. Keep the test-tube in a refrigerator until use.

However, as a general rule, it is better to subject all the specimens to microbiological analysis as soon as possible.

STORAGE

15-25 °C until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING and PRECAUTIONS

For laboratory use.

The product doesn't contain dangerous substances according to directives 1999/45/CE and 2001/60/CE or for which exist recognized exposure limits.

DISPOSAL of WASTE

Disposal of waste must be carried out according to national and local regulations in force.



Liofilchem s.r.l Bacteriology Products

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330

Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net



PRODUCT SPECIFICATIONS

NAME

STERILE SWABS WITHOUT ENVELOPE

PRESENTATION

Sterile dry swabs in tube.

STORAGE

15-25 °C.

PACKAGING

Code	Content	Packaging
85275	100 test	• 1 swab and 1 tube

USE

STERILE SWABS WITHOUT ENVELOPE are swabs for collection and transport of samples. They are constituted by sterile dry swabs inserted in a tube and endowed with a cap which hermetically closes the neck of the tube, so avoiding contamination of the tube and protecting user's hands. Sterilization is assured by irradiation at 2.5 Mrad rays.

TECHNIQUE

Refer to technical sheet of the product.

TABLE of SYMBOLS

Symbol	Meanings
	Catalogue number
	<i>In vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Kit content
	Use by
	Batch code
	Consult accompanying documents



Liofilchem s.r.l Bacteriology Products

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330
Web site: <http://www.liofilchem.net> E-mail: liofilchem@liofilchem.net