

Freeze-dried selective and growth supplements



Description	Packaging	Ref.
AMPICILLIN supplement (Ampicillin) Selective supplement for <i>Aeromonas hydrophila</i> isolation.	10 vials	81001
BACILLUS cereus supplement (Polymyxin B) Selective supplement for <i>Bacillus cereus</i> isolation (ISO 7932).	10 vials	81016
BORDETELLA supplement (Cephalexin) Selective supplement for <i>Bordetella</i> isolation.	10 vials	81013
BRUCELLA supplement (Polymyxin B, Bacitracin, Cycloheximide, Nalidixic acid, Nystatin, Vancomycin) Selective supplement for <i>Brucella</i> spp. isolation.	10 vials	81003
CAMPYLOBACTER growth supplement (Sodium pyruvate, Sodium metabisulphite, Ferrous sulphate) Growth supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81050
CAMPYLOBACTER Blaser Wang supplement (Polymyxin B, Trimethoprim, Amphotericin B, Cephalexin, Vancomycin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81051
CAMPYLOBACTER Butzler supplement (Bacitracin, Cycloheximide, Colistin sulphate, Sodium cephalozin, Novobiocin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81015
CAMPYLOBACTER C.C.D.A. supplement (Cefoperazone, Amphotericin B) Selective supplement for <i>Campylobacter</i> spp. isolation (ISO 10272-1).	10 vials	81037
CAMPYLOBACTER C.T.V.N. supplement (Cefoperazone, Trimethoprim, Vancomycin, Cycloheximide) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81038
CAMPYLOBACTER Karmali supplement (Sodium pyruvate, Cefoperazone, Cycloheximide, Vancomycin) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81036
CAMPYLOBACTER Preston supplement (polymyxin B, Rifampicin, Cycloheximide, Trimethoprim) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81004
CAMPYLOBACTER Skirrow supplement (Polymyxin B, Vancomycin, Trimethoprim) Selective supplement for <i>Campylobacter</i> spp. isolation.	10 vials	81055
C.F.C. (Pseudomonas) supplement (Cetramide, Fucidin, Cephalotin) Selective supplement for <i>Pseudomonas</i> isolation.	10 vials	81049
CHLORAMPHENICOL supplement (Chloramphenicol) Selective supplement for yeasts and moulds isolation.	10 vials	81017
CLOSTRIDIUM difficile supplement (D-Cycloserine, Cefoxitin) Selective supplement for <i>Clostridium difficile</i> isolation.	10 vials	81007
CLOSTRIDIUM perfringens (T.S.C.) supplement (D-Cycloserine) Selective supplement for <i>Clostridium perfringens</i> isolation (ISO 7937).	10 vials	81011
C.N. (Pseudomonas) supplement (Cetrimide, Nalidixic acid) Selective supplement for <i>Pseudomonas aeruginosa</i> isolation (ISO 12780).	10 vials	81006
C.N.A. (Staf / Strep) supplement (Colistin sulphate, Nalidixic acid) Selective supplement for staphylococci, streptococci, pneumococci isolation.	10 vials	81048
DERMATOPHYTE supplement (Chlortetracycline) Selective supplement for dermatophytes fungi isolation.	10 vials	81025
GARDNERELLA vaginalis supplement (Gentamicin, Nalidixic acid, Amphotericin B) Selective supplement for <i>Gardnerella vaginalis</i> isolation.	10 vials	81040
GENTAMICIN supplement (Gentamicin) Selective supplement for yeasts and moulds isolation.	10 vials	81033
HAEMOPHILUS supplement (Bacitracin, Vancomycin, Clindamycin) Selective supplement for <i>Haemophilus</i> spp. isolation.	10 vials	81014
KANAMYCIN supplement (Kanamycin) Selective supplement for group D streptococci isolation (APHA).	10 vials	81030

Freeze-dried selective and growth supplements



Description	Packaging	Ref.
KANAMYCIN / POLYMYXIN B supplement (Kanamycin, Polymyxin B) Selective supplement for <i>Clostridium perfringens</i> supplement.	10 vials	81031
L.C.A.T. supplement (Lincomycin, Colistin sulphate, Amphotericin B, Trimethoprim) Selective supplement for pathogenic neisseriae isolation.	10 vials	81012
LEGIONELLA growth supplement (L-Cysteine, Ferric pyrophosphate) Growth supplement for <i>Legionella</i> spp. isolation.	10 vials	80056
LEGIONELLA (B.C.Y.E.) growth supplement (L-Cysteine, Ferric pyrophosphate, ACES, Potassium hydroxide, Bromothymol blue, α-ketoglutarate) Growth supplement for <i>Legionella</i> spp. isolation (ISO 11731-1).	10 vials	81056
LEGIONELLA (B.M.P.A.) supplement (Anysomycin, Cefamandole, Polymyxin B) Selective supplement for <i>Legionella</i> spp. isolation.	10 vials	81002
LEGIONELLA (G.V.P.C.) supplement (Cycloheximide, Glycine, Polymyxin B, Vancomycin) Selective supplement for <i>Legionella</i> spp. isolation (ISO 11731-1).	10 vials	81008
LEGIONELLA (M.W.Y.) supplement (Anysomycin, Glycine, Polymyxin B, Vancomycin, Bromothymol blue, Bromocresol purple) Selective supplement for <i>Legionella</i> spp. isolation.	10 vials	81019
LISTERIA FRASER supplement (Ferric ammonium citrate) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81028
LISTERIA OXFORD supplement (Cefotetan, Fosfomycin) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81027
LISTERIA OXFORD WITH CEF/FOS/COL/CYC/ACR supplement (Cefotetan, Fosfomycin, Colistin, Cycloheximide, Acriflavine) Selective supplement for <i>Listeria</i> isolation (ISO 11290-1).	10 vials	81629
LISTERIA PALCAM supplement (Polymyxin B, Ceftazidime, Acriflavine HCl) Selective supplement for <i>Listeria monocytogenes</i> isolation (ISO 11290-1).	10 vials	81026
LISTERIA SUPPLEMENT FOR ENRICHMENT (FDA/IDF) Selective supplement for <i>Listeria</i> isolation.	10 vials	81045
MUG supplement (Methylumbelliferyl-β-D-glucosonide) Supplement for <i>Escherichia coli</i> identification.	10 vials	81020
O.N.P.G. 1.5% supplement (α-nitrophenyl-β-D-galactopyranoside) Supplement for enzyme β -galactosidase determination.	10 vials	81032
m-CP supplement Selective supplement for <i>Clostridium perfringens</i> isolation (D.L.vo 2/2/2001 np.31).	10 vials	81034
NOVOBIOCIN supplement (Novobiocin) Selective supplement for <i>Salmonella</i> isolation (ISO 6579).	10 vials	81021
ROSOLIC ACID supplement (Rosolic acid) Selective supplement for faecal coliforms isolation.	10 vials	81029
R.P.F. supplement (Rabbit plasma, Bovine fibrinogen, Tripsin inhibitor, Potassium tellurite) Selective supplement for coagulase positive staphylococci identification (ISO 6888-2).	8 vials	81057
SCHAEDLER supplement (Kanamycin, Vancomycin) Selective supplement for anaerobes isolation.	10 vials	81054
V.C.A.T. supplement (Vancomycin, Colistin sulphate, Amphotericin B, Trimethoprim) Selective supplement for pathogenic neisseriae isolation.	10 vials	81041
V.C.N. supplement (Vancomycin, Colistin sulphate, Nystatin) Selective supplement for pathogenic neisseriae isolation.	10 vials	81022
V.C.N.T. supplement (Vancomycin, Colistin sulphate, Nystatin, Trimethoprim) Selective supplement for pathogenic neisseriae isolation.	10 vials	81024
VITALEX growth supplement (Glutamine, Adenine, NAD) Growth supplement for fastidious microorganisms isolation.	10 vials	81023
YERSINIA supplement (Cefsulodin, Irgasan, Novobiocin) Selective supplement for <i>Yersinia</i> spp. isolation (ISO 10273).	10 vials	81039