



Self-contained for ethylene oxide sterilization		# Product with minimum request
Description	Packaging	Ref.
STERILtest BAT E6 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/disc).	100 units	91150
STERILtest BAT E5 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/disc).	100 units	91151 #