


Swabs



Description	Packaging	Ref.
STERILE DRYSWAB, FINE TIP, NARROW PLASTIC SHAFT Shaft: Flexible plastic for narrow orifices (urethral, pediatric). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85283
STERILE DRYSWAB, STRAIGHT ALUMINUM WIRE URETHRAL ENT Shaft: Controlled sampling from narrow orifices, ideal for male urethral sampling (urethral, ear). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85281
STERILE DRYSWAB, STANDARD PLASTIC SHAFT Shaft: standard plastic, suitable for most specimens (wound, nasal, throat, vaginal, skin). Alginate fiber: bud can be dissolved in weak organic acids to give total release of collected organisms, ideal for recovering organisms present in low numbers.	100 units	85282
STERILE DRYSWAB, ULTRA FINE FLEXIBLE TWIST-WIRE NASOPHARYNGEAL PEDIATRIC Shaft: ultra fine, twisted wire, for sampling nasopharyngeal region at back of nose (especially for whooping cough). Rayon fiber: non-inhibitory to fastidious microorganisms.	100 units	85280

Miscellaneous sterile disposable labware

Description	Packaging	ref.
Petri dishes 90 mm	600 units	96706
Petri dishes 60 mm	500 units	96705
Contact petri dishes 55 mm	500 units	96707
TWO SECTORS Petri dishes 90 mm	360 units	96709
THREE SECTORS Petri dishes 90 mm	500 units	96716
DISPOSABLE PLASTIC LOOPS 1 µL	100 units	96702
DISPOSABLE PLASTIC LOOPS 10 µL	100 units	96701
STERILE CONTAINER WITH SAMPLERS FOR STOOL COLLECTION	120 units	96708
TABLETS FOR CO ₂ GENERATION FOR ANAEROBOSIS	50 units	96740
DISPOSABLE PLASTIC PASTEUR PIPETTES	100 units	96703
TEST TUBES 16 mm x 160 mm (glass)	299 units	96717
CUPS FOR TEST TUBES (aluminum)	300 units	96718