

Tests for microbial identification		
Description	Packaging	Ref.
AESCULIN BILE TEST Rapid test for aesculin hydrolysis.	30 tubes	88008
ARGININE DECARBOXYLASE TEST Rapid test for arginine decarboxylation.	30 tubes	88016
BACITRACIN TEST Bacitracin 0.05 IU impregnated diagnostic discs for group A streptococci identification.	100 discs	9502
BETA LACTAMASE STICKS Sticks impregnated with Nitrocefin, used for the detection of β -Lactamase producing microorganisms.	30 sticks	88033
Brilliant Green 100 μg Discs for the presumptive identification of anaerobic bacteria.	50 discs	88041
CAMP TEST – R Rapid test for listeriae identification.	30 discs	88027
CAMP TEST – S Rapid test for listeriae identification.	30 discs	88021
CATALASE / OXY TEST Rapid test for catalase and peroxidase enzyme determination.	30 tubes	88023
H₂O₂ Reagent Reagent for catalase test.	5 mL 20 x 1 mL	80057 87003
CITRATE Disc for the identification of Enterobacteriaceae and non-fermenting gram-negative bacteria.	30 discs	88042
COAGULASE TEST Citratated freeze-dried rabbit plasma for coagulase test in tube.	5 vials (40 Tests)	88030
C-390 Discs containing antibiotic C-390 for the identification of <i>Pseudomonas aeruginosa</i> .	50 discs	88040
E. COLI TEST Rapid test for <i>Escherichia coli</i> identification (MUG, ONPG, Indole).	30 tubes 3 mL Kovac's	88006
GRAM TEST STICK Rapid test for gram-negative and gram-positive germs identification.	30 sticks 3 mL DMAC	88031
H₂S RAPID TEST Rapid test for H ₂ S production.	30 tubes	88013
HIPPURATE TEST Rapid test for hippurate hydrolysis.	30 tubes	88007
INDOLE TEST STICK Rapid test for indole production.	30 sticks 3 mL DMAC	88032
INDOLE TEST Rapid test for indole production.	30 tubes 3 mL Kovac's	88017
LISTERIA-mono TEST Rapid test for <i>Listeria monocytogenes</i> identification.	20 tubes 20 discs 3 mL reagent	88010
LYSINE DECARBOXYLASE TEST Rapid test for lysine decarboxylation.	30 tubes	88014
METRONIDAZOLE TEST Discs for <i>Gardnerella vaginalis</i> identification.	100 discs	9508
NINHYDRIN 7% Reagent for Hippurate Test.	10 mL bottle	80273
NITRATE TEST Rapid test for nitrates to nitrites reduction.	30 tubes 3 mL α -Naphthylamine 3 mL Sulphanilic acid	88009
O129 disc Test for differentiating <i>Vibro</i> spp. from other Gram negative bacteria.	50 discs	88043
ONPG Discs Rapid test for detecting the beta-galactosidase activity.	30 discs	88105 #
ONPG TEST Rapid test for lactose fermenting microorganisms differentiation.	30 tubes	88005
OPTOCHINE TEST Diagnostic discs for pneumococci identification.	100 discs	9501
ORNITHINE DECARBOXYLASE TEST Rapid test for ornithine decarboxylation.	30 tubes	88015
OXIDASE TEST DISCS Rapid test for cytochrome-oxydase enzyme detection.	30 discs	88004
OXIDASE TEST STICKS Rapid test for cytochrome-oxydase enzyme detection.	50 sticks	88029
OXIDASE TEST SWABS Rapid test for cytochrome-oxydase enzyme detection.	30 swabs	88003

Tests for microbial identification		
Description	Packaging	Ref.
PEPTIDASE A STICKS Rapid test for group A and group D streptococci identification.	30 sticks 5 mL Pyr reagent	88034
S.F. RAPID TEST Rapid test for group D streptococci determination.	30 tubes	88020
SULPHONAMIDE Diagnostic discs for <i>Gardnerella vaginalis</i> identification.	100 discs	9511
UREA / INDOLO TEST Rapid test for urea and indole detection.	30 tubes 3 mL Kovac's	88024
UREA RAPID TEST Rapid test for urease determination.	30 tubes	88011
V FACTOR TEST Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9504
X FACTOR TEST Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9503
X + V FACTOR TEST Diagnostic discs for <i>Haemophilus</i> spp. identification.	100 discs	9505